




Boyden Chamber Assay with CollOvine™

Collagen is a major component of the basement membrane and tissue scaffolding protein. The ability of tumor cells to invade through a collagen barrier is directly correlated with metastatic potential. The Boyden Chamber assay (**Figure 1**) is one of the most widely accepted cell migration study techniques.

Methods

Cell culture inserts with 8 μm pores in the PET membrane (ThinCert, Greiner Bio-One; 24-well size) were coated with 100 μL of CollOvine, the chambers were dried overnight at 37°C. HT-1080 cells (ATCC CCL-121, human fibrosarcoma cell line) or negative control cells (HEK-293) were deprived in serum-free culture medium the night before the migration experiment. On the day of the experiment, 700 μL of complete media (*i.e.*, with serum) was added to the bottom chamber, and 200 μL of cell suspension (6.8 x 10⁵ cells/mL) in serum-free media was added to the top chamber. The chambers were incubated overnight at 37°C and 5% CO₂. The next day, the culture media in the lower chamber were replaced with 450 μL serum-free media and 8 μM Calcein-AM and incubated for 45 minutes at 37°C and 5% CO₂. The inserts



were then transferred to fresh 24-well plates containing $500 \,\mu\text{L}\,0.2X$ Trypsin-EDTA per well and incubated for 10 minutes at 37°C and 5% CO₂, with agitation. Two hundred and fifty microliters of each detached cell suspension were then transferred to 96-well plates for quantification with a fluorescence plate reader at 485 nm excitation and 520 nm emission.




Figure 2. Fluorescence of HT-1080 cells migrated through a CollOvine-coated Boyden chamber.

Results

Figure 2 shows that the HT-1080 was able to invade through a CollOvine-coated membrane, while the negative control cells (HEK-293) was not able to invade through, as expected.