

CollOvine-PEG-MAL to rAspf9

Introduction

Maleimide conjugation chemistries are commonly used to generate covalent conjugates with sulfhydryl-containing molecules. The primary amine groups (lysine residues) of CollOvine are reacted with a heterobifunctional crosslinker to yield collagen derived with a maleimide group (CollOvine-PEG2-MAL), which is then reactive with other species containing free sulfhydryl groups (*e.g.*, cysteine residues). Cysteine residues commonly occur naturally in protein sequences or can easily be engineered into the protein or peptide sequence.

We have successfully conjugated CollOvine-PEG2-MAL to rAspf9, an *Aspergillus fumigatus* cell wall glucanase protein. The rAspf9 was recombinantly prepared in-house in an *E. coli* expression vector. The resulting protein is ~29 kDa, with one cysteine residue close to the N-terminus. The rAspf9 protein has been shown to be an effective antigen in a developmental *A. fumigatus* vaccine (AspaVax) by our collaborative partner, Molecular Express, Inc.

Procedure

To prepare the cysteines of the rAspf9 for conjugation, one volume of rAspf9 at ~ 2 mg/mL was reduced with 0.1 volume of 2-mercaptoethanol for ~ 30 minutes at room temperature. The excess 2-mercaptoethanol was removed by passage through a desalting column (e.g., Illustra NAP-5 column). The peak fractions




Figure 1. SDS-PAGE and Western Blot of CollOvine-rAspf9 conjugation. Lane 1: Protein Marker; Lane 2: CollOvine-PEG2-MAL and rAspf9; Lane 3: CollOvine-PEG2-MAL only; Lane 4: rAspf9 only. The image on the left is the Coomassie-stained SDS-PAGE analysis of the samples indicated. The image of the right is the Western Blot against rAspf9 of the samples indicated.

were determined by absorbance at 280 nm. The reduction of the cysteine residues was confirmed by Ellman's reagent. To avoid re-oxidation of the cysteines, the reduced rAspf9 protein was used immediately.

One mole of CollOvine-PEG2-MAL was mixed with about 5 mole excess of reduced rAspf9 protein. The reaction was incubated for 30 minutes at room temperature. The conjugates were analyzed as described below.

Results

The conjugate was analyzed by SDS-PAGE (**Figure 1, left**) and Western Blot (**Figure 1, right**). CollOvine-PEG2-MAL, rAspf9 and CollOvine-PEG2-MAL reacted with rAspf9 were loaded on an SDS-PAGE gel. One gel was stained with Coomassie stain to indicate the molecular weights of the samples (**Figure 1, left**). Another gel with the same samples was analyzed by Western Blot (**Figure 1, right**) following a standard procedure. The primary antibodies were anti-rAspf9 antibodies obtained from the serum of mice vaccinated with AspaVax, the secondary antibodies were goat anti-mouse IgG1 antibodies conjugated to horse radish peroxidase (ThermoFisher Scientific), and the substrate was the chromogenic CN/DAB substrate (ThermoFisher Scientific).



Technical Bulletin

Updated on January 12, 2021

As shown in lane 2 of the SDS-PAGE (**Figure 1**, **left**), the CollOvine-PEG2-MAL-rAspf9 conjugate has a higher molecular weight than the CollOvine-PEG2-MAL alone (lane 3) or rAspf9 alone (lane 4). Since an excess of rAspf9 was added, unconjugated rAspf9 (~30 kDa band) remained in lane 2. The Western Blot results (**Figure 1**, **right**) show that the anti-rAspf9 antibodies are reactive against the rAspf9 (see lanes 2 and 4, ~30 kDa). The dark color at the higher molecular weight in lane 2 (**Figure 1**, **right**) confirms the presence of CollOvine-PEG2-MAL-rAspf9, and that the CollOvine-PEG2-MAL was conjugated with multiple copies of the rAspf9 protein.