

Gelation Procedure for CollOvine-RD (3 mg/mL) Lot Number: 190919005FS010

Materials:

- 1. CollOvine-RD (3 mg/mL, Lot 190919005FS010)
- 2. 1 M TRIS Buffer, pH 8

Procedure:

- 1. To an appropriately sized container, transfer the amount of CollOvine-RD needed.
- 2. Add 80 µL of 1 M TRIS Buffer to every 1 mL of CollOvine-RD.
- 3. Mix the solution **thoroughly**, but gently to avoid bubbles. Do not vortex.
- 4. If needed, aliquot the mixture into desired container(s).
- 5. Incubate the neutralized mixture at 37°C to initiate gelation. Allow the gel to form for at least 3 hours. Avoid moving the mixture during the incubation time.