

CollOvine-RDP Rehydration Protocol

Materials Provided:

- 1. 600 mg CollOvine-RDP
- 2. 125 mL of 0.02 M acetic acid (sterile), store at room temperature

Materials needed but not provided

- 1. Sterile pipet (25 or 50 mL) or equivalent
- 2. Pipet aid
- 3. Orbital shaker or equivalent
- 4. Sterile room or biosafety cabinet (laminar flow hood)

Procedure:

The lyophilized collagen must be hydrated with small incremental amounts of 0.02 M acetic acid. OviGenex provided 0.02 M acetic acid for your convenience. The following procedure is for the hydration of the bulk lyophilized powder to 6 mg/mL final concentration.

- 1. To maintain the sterility of the collagen, only open the bottle under sterile conditions (*i.e.*, a biosafety cabinet).
- 2. To the 600 mg of lyophilized CollOvine-RDP, add 50 mL of 0.02 M acetic acid. Swirl the bottle to mix and place on an orbital shaker for 4 hours at room temperature.
- 3. Add 25 mL of 0.02 M acetic acid and continue mixing on an orbital shaker for 2 hours at room temperature. Store solution at 2 8 °C overnight.
- 4. Following day, add 24 mL of 0.02 M acetic acid and continue mixing on an orbital shaker for 8 hours at room temperature until dissolved.
- 5. Store the collagen solution at 2 8 °C.