



## PRODUCT DESCRIPTION SHEET

CollOvine-RD is purified atelocollagen from the dermis of prion-free sheep raised in Western Australia.

**Product name:** CollOvine-RD

**Product code:** 1040116700

**Lot no:** 190919005FS010

<b>Characterization Method</b>	<b>Results</b>
Appearance (Form)	Translucent liquid
Collagen Concentration in 0.02 M acetic acid	3.0 mg/mL
Purity (SDS-PAGE)	≥ 99 %
Density (g/mL)	1.01 (+/- 0.02)
pH	3.5 (+/- 0.25)
Total aerobes	None detected
Yeast	None detected
Mold	None detected
Coliform	None detected
Collagen Identity with Trypsin resistance	Characteristic banding pattern on SDS-PAGE
Endotoxin (EU/g collagen)	≤ 100 EU/g