

PRODUCT DESCRIPTION SHEET

CollOvine-RDP is purified at elocollagen from the dermis of prion-free sheep raised in Western Australia.

Product name: CollOvine-RDP **Product code**: 1041116703 **Lot no**: 190919005L001

Characterization Method	Results
Appearance (Form)	White powder cake
Collagen	600 mg collagen
Purity (SDS-PAGE)	≥ 99 %
Total aerobes Yeast Mold Coliform	None detected None detected None detected
Collagen Identity with Trypsin resistance	Characteristic banding pattern on SDS-PAGE
Endotoxin (EU/g collagen)	$\leq 100 \text{ EU/g}$

Phone: 562-297-0310 / Fax: 310-635-5503